



Deltaclon

For research use only

COVID-19 SARS-CoV-2 2019-nCoV

| Proveedor | Cat. | Descripción |
|------------------|-----------------------------|---|
| Absource | CD302-02 | 2019-nCoV Triplex RT-qPCR Detection Kit 100 test |
| Absource | C6603C | 2019-nCoV IgG/IgM Detection Kit (Colloidal Gold-Based) 50 test |
| AccuBiotech | AB-IDT-B352 | Accu-Tell COVID-19 IgG/IgM Cassette (Whole Blood/Serum/Plasma) 20T/40T |
| Acro Biosystems | AC2-H82E6 | Human HEK293 Biotinylated Human ACE2 / ACEH Protein, His,Avitag™ |
| Acro Biosystems | SPN-C52H4 | 2019-nCoV HEK293 2019-nCoV S protein Full Length (R683A, R685A), His Tag |
| Acro Biosystems | S1N-C52H3 | 2019-nCoV HEK293 2019-nCoV S1 protein, His Tag |
| Alpha Diagnostic | RV-404100-1 | Human Anti-SARS-CoV-2 Virus (COVID-19) Spike 1 IgG ELISA Kit |
| Biologix | | Virus Transport and Preservation Medium |
| Biologix | | COVID-19 Nucleic acid detection kit (multiplex PCR fluorescent probe method), |



Deltaclon SL

C/ Javier del Quinto, nº 4 bajo 3

28043 Madrid

ES-B81380370

Teléfono 91 764 53 64

Fax 91 381 42 75

info@deltaclon.com

www.deltaclon.com



Deltaclon

For research use only

| | | |
|----------|-------|---|
| Biologix | | A Rapid IgM-IgG Combined Antibody Test Kit |
| Biologix | W195 | SARS-CoV-2 Antibody Test (Lateral Flow Method) |
| Biotnt | | Realtime RT-PCR COVID-19 KIT |
| Biovison | A2069 | 100 µg Anti-Canine Coronavirus (Capture) Antibody |
| Biovison | A2068 | 100 µg Anti-Canine Coronavirus (Self-pair) Antibody |
| Biovison | A2070 | 100 µg Anti-Canine Coronavirus Polyclonal Antibody |
| Biovison | A2063 | 100 µg Anti-MERS & SARS-CoV NP Mouse IgG1 Ab |
| Biovison | A2065 | 100 µg Anti-MERS-CoV S1 Mouse IgG1 Antibody |
| Biovison | A2067 | 100 µg Anti-MERS-CoV S1 Mouse IgG2a Antibody |
| Biovison | A2066 | 100 µg Anti-SARS-CoV NP Mouse IgG1 Antibody |
| Biovison | A2064 | 100 µg Anti-SARS-CoV NP Mouse IgG2b Antibody |
| Biovison | A2061 | 50 µg Anti-SARS-CoV-2 Antibody |
| Biovison | A2060 | 50 µg Anti-SARS-CoV-2 Antibody (Clone# 6F10) |
| Biovison | B2997 | 500 µg, 1 mg Anti-viral Agent Remdesivir |
| Biovison | K1460 | 100 Rxns Kit Coronavirus (SARS-CoV-2) PCR Detection Kit |
| Biovison | P1503 | 10 µg, 50 µg Envelope Protein (SARS-CoV ENV;1-76) |
| Biovison | P1504 | 10 µg, 50 µg Matrix Protein (SARS-CoV; 182-216) |
| Biovison | P1505 | 10 µg Membrane Protein (SARS-CoV, His tag) |

Deltaclon SL

C/ Javier del Quinto, nº 4 bajo 3

28043 Madrid

ES-B81380370

Teléfono 91 764 53 64

Fax 91 381 42 75

info@deltaclon.com

www.deltaclon.com



Deltaclon

For research use only

| | | |
|----------|------------|---|
| Biovison | P1506 | 10 µg, 50 µg Nucleoprotein (CoV-NP 229E) |
| Biovison | P1507 | 10 µg, 50 µg Nucleoprotein (CoV-NP-NL63) |
| Biovison | P1509 | 10 µg Nucleoprotein (SARS-CoV) |
| Biovison | P1510 | 10 µg, 50 µg Nucleoprotein (SARS-CoV; 1-49) |
| Biovison | P1511 | 10 µg, 50 µg Nucleoprotein (SARS-CoV; 1-49, 192-220) |
| Biovison | P1508 | 10 µg, 50 µg Nucleoprotein (SARS-CoV; 340-390) |
| Biovison | P1512 | 10 µg Nucleoprotein (SARS-CoV; 340-390) |
| Biovison | P1513 | 10 µg, 50 µg Spike Protein (MERS-CoV S1; 18-725) |
| Biovison | P1514 | 10 µg, 50 µg Spike Protein (MERS-CoV S1; 56-295) |
| Biovison | P1520 | 10 µg Spike Protein (SARS-CoV S, His tag) |
| Biovison | P1515 | 10 µg, 50 µg Spike Protein (SARS-CoV S1; 12-53, 90-115, 171-203) |
| Biovison | P1516 | 10 µg Spike Protein (SARS-CoV S1; His tag) |
| Biovison | P1519 | 10 µg, 50 µg Spike Protein (SARS-CoV S2) |
| Biovison | P1518 | 10 µg, 50 µg Spike Protein (SARS-CoV S2; 408-470, 540-573) |
| Biovison | P1517 | 10 µg Spike Protein (SARS-CoV S2; 408-470, 540-573, His Tag) |
| Creative | Spike-191V | Recombinant COVID-19 Spike protein S1 subunit, His-tagged |
| Creative | Spike-195V | Recombinant COVID-19 Spike protein S1 subunit, His-tagged |
| Creative | DAGC089 | Recombinant SARS-CoV-2 Spike Protein Receptor Binding Domain [Fc] ELISA, LFIA, WB Antigen |
| Creative | DAGC091 | Recombinant SARS-CoV-2 Spike Protein 1 [His] SDS-PAGE Antigen |



Deltaclon SL

C/ Javier del Quinto, n° 4 bajo 3

28043 Madrid

ES-B81380370

Teléfono 91 764 53 64

Fax 91 381 42 75

info@deltaclon.com

www.deltaclon.com



Deltaclon

For research use only

| | | |
|-----------|------------|---|
| Creative | DAGC092 | Recombinant SARS-CoV-2 Spike Protein 1 [human Fc] SDS-PAGE Antigen |
| Creative | DAGC093 | Recombinant SARS-CoV-2 Spike Protein 1 [mouse Fc] SDS-PAGE Antigen |
| Creative | DAGC094 | Recombinant SARS-CoV-2 Nucleocapsid Protein [His] ELISA, LFIA, WB Antigen |
| Creative | DAGS029 | Recombinant SARS-CoV-2 Nucleocapsid Protein (29-213 aa) [His] ELISA, LFIA, WB Antigen |
| Creative | DAGS059 | Recombinant SARS-CoV-2 Nucleocapsid Protein (59-213 aa) [His] ELISA, LFIA, WB Antigen |
| Creative | CABT-RM320 | Rabbit anti-SARS-CoV-2 NP monoclonal antibody, clone 120 ELISA, WB Antibody |
| Creative | CABT-RM321 | Rabbit anti-SARS-CoV-2 spike glycoprotein monoclonal antibody, clone 118 ELISA Antibody |
| Creative | CABT-CS026 | Mouse Anti-SARS-CoV-2 NP Monoclonal antibody, clone 7G21 ELISA, WB, IHC Antibody |
| Creative | CABT-CS027 | Mouse Anti-SARS-CoV-2 NP Monoclonal antibody, clone 4B21 ELISA, IHC Antibody |
| Creative | CABT-CS024 | Rabbit Anti-SARS-CoV-2 NP Polyclonal antibody ELISA Antibody |
| Creative | DEIA1035 | SARS Coronavirus IgG ELISA Kit |
| Creative | DTS920 | 2019-nCoV IgG/IgM Rapid Test (cassette) 200 T (MOQ) |
| Creative | DEIA2021 | COVID-19 Nucleoprotein ELISA Kit |
| Creative | DEIA2022 | COVID-19 Spike glycoprotein ELISA Kit |
| Creative | CD019RT | 2019-Novel Coronavirus Multiplex RT-qPCR Kit duplex (two pairs of PCR primers) 50 T |
| Creative | CD019RT | 2019-Novel Coronavirus Multiplex RT-qPCR Kit triplex (three pairs of PCR primers) 50 T |
| Creative | DEIA2023 | SARS-CoV-2 IgG ELISA Kit |
| Dynamiker | DNK-1419-1 | 2019-nCoV IgG/IgM Rapid Test (LFA) |
| Dynamiker | DNK-1418-1 | Novel Coronavirus (2019-nCoV) RT-PCR Kit |

Deltaclon SL

C/ Javier del Quinto, n° 4 bajo 3

28043 Madrid

ES-B81380370

Teléfono 91 764 53 64

Fax 91 381 42 75

info@deltaclon.com

www.deltaclon.com



Deltaclon

For research use only

| | | |
|--------------------|-----------|---|
| Eosmedchem | | Remdesivir GS-5734, CAS 1809249-37-3 |
| Eosmedchem | | COVID-19 Rapid Test Cassette, CE certification |
| Epitopediagnostics | | EDI™ Novel Coronavirus COVID-19 IgG ELISA Kit |
| Epitopediagnostics | | EDI™ Novel Coronavirus COVID-19 IgM ELISA Kit |
| GeneScrip | L00831 | SARS-CoV-2 Spike S1-RBD IgG&IgM ELISA Detection Kit |
| Glory Science | | Coronavirus COVID-19 Antibody Detection Kit |
| Headwaychina | | Corona virus rapid test kits, with CE certificate. |
| Honorbiotech | | 1 COVID-19 IgG/IgM Rapid Lateral Flow Test (colloidal gold-based) |
| Honorbiotech | | COVID-19 ORF1ab Gene Isothermal Nucleic Acid Amplification Kit (RT-SRDA Fluorescent Type) |
| Honorbiotech | | COVID-19 N Gene Isothermal Nucleic Acid Amplification Kit (RT-SRDA Fluorescent Type) |
| Honorbiotech | | COVID-19 ORF1ab / N Dual Nucleic Acid Detection Kit (RT-qPCR method). |
| Honorbiotech | | 3 COVID-19 (SARS-CoV-2) Nucleocapsid Protein |
| Honorbiotech | | 4 COVID-19 (SARS-CoV-2) Spike (full S sequence) Recombinant antigen |
| Honorbiotech | | 5 COVID-19 (SARS-CoV-2) RDB Recombinant Antigen |
| Sunlong | SL-COV | COVID-19 IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma) |
| TRC | R530520 | Remdesivir Carboxylic Acid |
| Vazyme | C6603C-02 | 2019-Novel Coronavirus (2019-nCoV) IgG / IgM Detection Kit (Colloidal Gold-Based) 50 test |
| Vazyme | CD302-02 | 2019-Novel Coronavirus (2019-nCoV) Triplex RT-qPCR Detection Kit |



Deltaclon SL

C/ Javier del Quinto, n° 4 bajo 3

28043 Madrid

ES-B81380370

Teléfono 91 764 53 64

Fax 91 381 42 75

info@deltaclon.com

www.deltaclon.com